



ZEPF Augmentation Material Applicator

19.714.16 - ZEPF Augmentation Material Applicator Ø 5mm straight, with a lateral opening to insert the augmentation materia



41.868.07 - Kirsch Sinus 7 spoon instrument ZEPF-Titanium-Line, double-ended, blue, 17,5cm, spoon 8x10 mm



ZEPF Augmentation Material Applicator

The ZEPF augmentation material applicator supports the fractional filling of the subantral space with particulate augmentation material. To compact the augmentation material you use a sinus plugger.



The instrument with diameter of 5mm is loaded laterally.



The augmentation material can be applied precisely with the injecting mechanism.

The instrument can be taken apart for optimal cleaning.

NEWS

03/09



HELMUT ZEPF
MEDIZINTECHNIK GMBH